

 enova®



**SEE BETTER  
PERFORM  
BETTER**  
ZUMAX® Microscopes

**zumax**

# REVOLUTIONIZE YOUR PRACTICE WITH A ZUMAX® MICROSCOPE

Magnification using the latest in microscope technology gives professionals like you the opportunity to offer a higher standard of care in dentistry. Enova Illumination is proud to be the exclusive distributor of Zumax® microscopes in the U.S. and Canada, bringing you the highest quality and the latest innovations.



## ENHANCED MAGNIFICATION

Increased precision in all aspects of patient care - from diagnosis and treatment to efficiency and patient outcomes.



## SUPERIOR ERGONOMICS

Reduce neck and back injuries that can limit a career by not relying on loupes and insufficient illumination. Improve your posture by utilizing heads-up magnification with a monitor. 3D visualization is optional.



## HIGH QUALITY OPTICS

On all Zumax® microscopes, all lenses are multi-coated, anti-reflective, apochromatic, and made from hand-selected German glass.



## INTEGRATED 4K CAMERA

Connect to a monitor or your network for the ability to stream and save high-definition photos and videos. Edit videos immediately for training and sharing with consulting physicians. Display real-time imaging for increased patient education.



## DOCUMENTATION AND TRAINING

Microscopes allow you to easily show, discuss and document treatments and procedures, both with patients and other specialists.



## RETURN ON INVESTMENT

Elevate your practice with Zumax® Microscopes, empowering you to confidently perform complex procedures with enhanced visual precision!

Scan here to learn more!



# Zumax® Dental Microscopes

## 3200 R2/PRO

### FEATURES

- **Joystick-Controlled Continuous Zoom:** Magnification ranging from ~1-20x AutoFocus is standard
- **Electromagnetic Clutch:** One-click magnetic clutch for ease of positioning and smooth movement
- **Integrated 4K Camera:** Seamlessly save, record and stream images/videos to your network and monitor(s)  
\*3D 4K camera optional
- **Integrated LED Illumination:** ≥105,000 lux at 250mm working distance, multiple spot sizes and light and laser filters included. LED lifetime: 60,000 hours
- **Integrated VarioScope:** Working range 200mm to 450mm
- **Mounting Options:** Floor Stand, Ceiling Mount, Wall Mount



**Motorized Continuous  
Zoom System**



**Electromagnetic  
Clutch**



**Integrated 4K  
Camera**

# Zumax<sup>®</sup> Dental Microscopes

## 2380 R2

### FEATURES

- **Continuous Zoom Magnification:**  
Ranges from ~1-20x
- **FreeFlow balancing arm:**  
Free flowing, ultra smooth movement to simplify positioning
- **Integrated 4K Camera:**  
Seamlessly save, record and stream images/videos to your network and monitor(s)  
\*3D 4K camera optional
- **Integrated LED Illumination:**  
≥100,000 LUX at 250mm working distance, multiple spot sizes and light and laser filters included. LED lifetime: 60,000 hours
- **Integrated VarioScope:**  
Variable working distance 200mm-450mm
- **Mounting Options:**  
Floor Stand, Ceiling Mount, Wall Mount



**Continuous  
Zoom System**



**VarioScope Fine  
Focus**



**Integrated 4K  
Camera**

# Zumax® Dental Microscopes

## 2050

### FEATURES

- **Continuous Zoom Magnification:**  
Ranges from ~1-20x
- **ErgoPendulum Arm:**  
Binoculars remain stationary as microscope pivots right and left
- **Integrated LED Illumination:**  
≥75,000 LUX at 250mm working distance, multiple spot sizes and light and laser filters included. LED lifetime: 60,000 hours
- **4K Camera Module (Optional):**  
Seamlessly save and display images and video
- **Mounting Options:**  
Floor Stand, Ceiling Mount, Wall Mount, Tabletop



**Continuous  
Zoom System**



**Integrated LED**

# Zumax® Dental Microscopes

## 2350

### FEATURES

- **6-Step Magnification:**  
Ranges from ~1-24x
- **Integrated LED Illumination:**  
≥65,000 LUX at 250mm working distance, multiple spot sizes and light and laser filters included. LED lifetime: 60,000 hours
- **Integrated HD Camera (Optional):** Seamlessly connects to your network
- **Integrated VarioScope:**  
Variable working distance 190mm-300mm
- **Mounting Options:**  
Floor Stand, Ceiling Mount, Wall Mount, Tabletop



**Integrated LED**



**VarioScope Fine Focus**

# Microscope Accessories

## BROKEN INSTRUMENT REMOVAL KIT

Instrument separation is one of the most common complications in the root canal treatment process. The innovative design of the Zumax® Broken Instrument Removal Kit makes instrument removal quick and easy.


Kit Includes:

- 3 Micro-Trephine Burs
- 1 Extractor Handle
- 3 Extractors
- 1 Foreign Forceps

Extractor Handle and Tip



Extractor Tip Options




## MICRO DENTAL TWEEZER

The Zumax® Micro Dental Tweezer is designed for accurate foreign material removal - no matter the angle or complex accessibility.

Kit Includes:

- 1 360 Tweezer Head
- 1 Handle
- 1 Cleaning Adapter



## MOBILE ADAPTER

Connect mobile devices to Zumax® microscopes to allow you to instantly share your images through live-streaming, text messages, or email.

Compatible with:  
Most smartphones via universal adapter  
iPad mini (6th Gen.)



# ENOVA SATISFACTION GUARANTEE

Use our products and love them. If not, we'll take them back. It's that simple.  
For full warranty specifications visit [www.goenova.com/product-support](http://www.goenova.com/product-support)

## Connect with us!



To learn more about how Enova can help you see better, so you can perform better, visit [www.goenova.com](http://www.goenova.com), call 651.236.8857, or email [sales@goenova.com](mailto:sales@goenova.com)

### World Headquarters

Enova Illumination  
6305 Sandburg Road, STE 100  
Minneapolis, MN 55427

ISO 13485: 2016 Certified QMS  
MDR Class 1 Certified

F20578C

© 2024 Enova Illumination, Inc. All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.  
Oakley® and Flak® are trademarks of Oakley, Inc.  
Zumax® is a registered trademark of Zumax Medical Co., LTD

