


M1 Bright



- Integrated LED
- 65,000 LUX Illumination
- Mixed Colors



M2 Rainbow



- Integrated Endo-Gauge
- Mixed Color Assortment



M2 Mag



- Integrated Magnification
- Integrated Endo-Gauge
- Available in 2x or 3x



M3 Breeze




- Integrated AirClean Anti-Fog
- Compatible with
 - Dental Unit
 - 3 Way Syringe
- Mixed Color Assortment



M5 Mini



- Integrated AirClean Anti-Fog
- Ø15mm Mirror
- Compatible with
 - Dental Unit
 - 3 Way Syringe
- Mixed Color Assortment



M7 Eva Suction



- Integrated Suction
- Integrated Suction
 - 21 Holes Provide Suction in All Directions
- Mixed Color Assortment





**MOTRANSER
DENTAL MIRRORS**

enova[®]
goenova.com/dental-mirrors

All Motranser Dental Mirrors Feature:

- **Ultra High Definition Mirrors: 98% Reflectance Rate**
- **Alloy Coated Mirrors: Scratch and Chemical Resistant**
- **Autoclavable Mirror Handles: 135° C for up to 500 cycles**
- **Color Coding for Easy Identification**

To learn more or to purchase your new mirrors, visit:

goenova.com/dental-mirrors

To learn more about how Enova can help you see better, so you can perform better, visit www.goenova.com, call 651.236.8857, or email sales@goenova.com.

World Headquarters

Enova Illumination
6305 Sandburg Road, STE 100
Minneapolis, MN 55427

F20634A
© 2025 Enova Illumination, Inc. All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.
Motranser® is a registered trademark of ACME TOOLING TECHNOLOGY CO., LTD.

enova[®]