

# QUBIT™ WARM LED HEADLIGHT



All Comfort Pads  
Treated with MICROBAN®  
Antimicrobial Technology 

Adjustable Headband  
for Customizable Fit

Adjustable Links  
for Ergonomic Wear

Durable  
Aluminum Housing

Fixed  
Spot Size

**Beam Adjustment Filter (optional)**

Doubles the otherwise fixed spotlight size.

## THE MOST ECONOMICAL LED HEADLIGHT

With a color-rendering index (CRI) of 95, the **Qubit Warm LED Headlight** features a warm light that offers the truest color rendition on the market. This headlight offers a lightweight and economical way to bring LED light into any practice. Attach the Beam Adjustment Filter to make the headlight semi-adjustable.

### 80,000 - 160,000 LUX

Lux ranging from 80,000 - 160,000 lux\* — depending on desired working distance — of pure homogenous light for true tissue color rendition.

### COMFORTABLE AND ERGONOMIC

Designed for ultimate adjustability, comfort, and durability, the Quantum Cool headlight adjusts to your exact fit, desired working distance and light angles.

### BEAM ADJUSTMENT

The simple flip filter allows the Qubit to have a semi-adjustable light that doubles the spot size. For smaller, more focused light, simply adjust the filter off the lens.

\*Lux is variable based on preferred working distance

# TECHNICAL DATA

Model Name	Quest	Quantum Warm	Quasar Warm	Qubit Warm	Qubit Pro
Light Source	High-powered LED	LED refractive technology			
Brightness @ 10"	240,000 Lux	350,000 Lux	230,000 Lux	160,000 Lux	120,000 Lux
Brightness @ 14"	125,000 Lux	175,000 Lux	115,000 Lux	80,000 Lux	55,000 Lux
Color Temp.	6,200K	4,500K	4,500K	4,500K	6,500K
Spot Size @ 14"	2.8-4.8 in	2-5 in	2-5 in	2 in Fixed	3.6 in Fixed
Headlight Housing	Thermoplastic, Aluminum	Aluminum	Aluminum	Aluminum	Aluminum
Headlight Weight	10.7oz/303g	12.3oz/348g	6.6oz/188g	6.6oz/186g	6.1oz/174g
Compatible Battery	PP4+ (Plus)	PP4		PP2 or PP4	
Battery Run-Time	PP4+ (Plus): 4-18 hrs	PP4: 5-10 hrs	PP2: 5-10 hrs, PP4: 10-24 hrs		
Battery Weight	PP4+ (Plus): 10.3oz/292g	PP4: 10.6oz/300g	PP2: 6.5oz/184g, PP4: 10.6oz/300g		
Battery Size	PP4+ (Plus): 3.2×4×1 in	PP4: 3.2×4×1 in	PP2: 2.3×4×1 in, PP4: 3.2×4×1 in		

# MODEL INFORMATION

	Quest	Quantum Warm	Quasar Warm	Qubit Warm	Qubit Pro
<b>FEATURES</b>					
Adjustable spot size	●	●	●	–	–
Pivot bracket	●	●	●	●	●
Hybrid thermal management using fan	●	–	–	–	–
Upgrade battery	–	–	●	●	●
Optional UV filter	–	–	●	●	●
<b>PACKAGE: DEFAULT vs OPTION</b>					
Single ratchet	●	●	●	●	●
Single ratchet w/occipital	●	●	●	●	●
Double ratchet	●	●	●	●	●
Double ratchet w/occipital	●	●	●	●	●
PP2 battery (Rechargeable Li-Ion)	–	–	●	●	●
PP4 battery (Rechargeable Li-Ion)	–	●	●	●	●
PP4 Plus battery (Rechargeable Li-Ion)	●	–	–	–	–
<b>WARRANTY</b>					
Headlamp covering manufacturer defects - 5 years	–	●	●	●	●
Headlamp covering manufacturer defects - 1 year	●	–	–	–	–
Headlamp pigtail - 1 year	●	●	●	●	●
Batteries and accessories - 1 year	●	●	●	●	●



Test the brilliance of the Enova headlight through our complimentary evaluation program. Connect with your sales representative now to experience surgery in a whole new light.

## Enova Illumination

O: 651.236.8857

E: sales@goenova.com

W: goenova.com

A: 6305 Sandburg Rd, STE 100, Minneapolis, MN 55427



F20438 Rev. E

Copyright © 2025 Enova Illumination, Inc. All rights reserved.

All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.

Microban® is a registered trademark of Microban Products Company.

Antimicrobial Technology does not protect the user against disease-causing pathogens.