

# proCAM XS BY FUTUDENT



Versatile mount is compatible with most loupes, headlights, and overhead lights

Optional customized light adapter

Adjust focus to your exact working distance



Lightweight: At just 14 grams, you'll barely notice a difference when mounted to your loupe or headlight

Durable Aluminum Housing

## LIGHTS, CAMERA, DOCUMENTATION

Modern healthcare presents an increasing need for safe, non-invasive and efficient video and photo capture for clinical documentation, lecturing, and education. Weighing just 14 grams, the **Futudent proCam XS** is a miniature camera with giant 4K features and easy mounting to loupes, headlights and overhead lights.

### FAST, SIMPLE EDITING

Proprietary Futudent software captures and instantly makes your video and photos available so you can focus on your patient, not your camera.

### HANDS-FREE OPERATION

Start and stop video recording and capture photos without lifting a finger. The proCam XS is powered by a footpedal so you don't have to interrupt treatment to document it.

### TREATMENT ACCEPTANCE

Increase your ROI and treatment acceptance rates by being able to instantly help patients visualize their case and comprehend your treatment plan.

# TECHNICAL DATA

**Resolution:** 4K30/1080p30/720p30/8Mpix stills

**Weight:** 14 grams (0,49 ounces)

**Dimensions:** 1.18" x 0.79" x 0.79" (30 × 20 × 20 mm)

**Optics:** Exchangeable UHD lenses with manual focus (8/16/25/35 mm)

**Construction:** Full Aluminum Body

**Connectivity:** USB 2 with detachable cable (USB micro connector)

**Video Editing Software:** Futudent software; 5 licenses included/unit

**Camera Warranty:** 2 years

**Compliance:** MDR 2017/745, Classification 1

| Mounting Options |                                                 |
|------------------|-------------------------------------------------|
| Loupe Mount      | Universal and Specific Loupe Adapters           |
| LED Headlight    | Headband and Headlight Adapters                 |
| Free Standing    | Overhead Light, Gooseneck/Flexible Arm Adapters |

## PAIRS WITH YOUR FAVORITE ENOVA PRODUCTS

The latest cameras from Futudent are so lightweight, you'll hardly notice the difference when you mount the camera on your loupe or headlight.



**Vizix® Galilean Loupes**



**Vizix® Prismatic Loupes**



**Quasar™ LED Headlight**



**Qubit™ LED Headlight**

To learn more about how Enova can help you see better, so you can perform better, visit [www.goenova.com](http://www.goenova.com), call 651.236.8857, or email [sales@goenova.com](mailto:sales@goenova.com).

### World Headquarters

Enova Illumination  
6305 Sandburg Road STE 100  
Minneapolis, MN 55427

F20474, Rev B  
© 2022 Enova Illumination, Inc. All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.

