

# ERGO ELITE

## VIZIX® LOUPES



Prescription-friendly telescopes and carrier lens to create a precision fit unique to your eyes

Every frame comes with a magnetic adapter to securely hold a light or camera

All-glass telescopes lenses for high quality optical resolution

Built with aluminum housings for easy cleaning and greater durability

## LIGHTWEIGHT AND ERGONOMIC:

The Enova **Vizix® Ergo Elite** Prismatic Loupe is the ultimate accessory to enhance your practice. Modern engineering, found in Enova Vizix® Loupes, places the telescope lower on the carrier lens for a greater range of vision. When used properly, in the right setting, it will help decrease fatigue and pain associated with repetitive cervical spine movement. The Ergo Elite is available in 3.5x, 4.5x and 5.5x! All carrier frames have been hand-selected by experts for ergonomics, comfort, and durability. Select from Enova's Airon® family of frames for a stylish presentation.

### GREATER RANGE OF VISION

Strategically placed telescopes provide an expansive range of vision, allowing you to easily work around equipment without having to take off your loupes.

### LIGHTER & MORE COMPACT

Telescope is designed to be lightweight, without sacrificing technology that supports decreasing stress on your neck, shoulders, and lower back for long-term body health.

### FUNCTIONAL FRAMES

Enova frames provide the complete package: modern style and superior durability. Comfort, ergonomics, and functionality are all seamlessly integrated into an attractive design.

# OCULAR SPECIFICATIONS

\*Oculars weighed as a pair without frame

	Galilean	Prismatic	Ergo Elite	Micro Ergo Elite
Image				
(Ocular Color)	 *3.5x is only available in Black/Silver			
Magnification	2.5x / 3.0x / 3.5x	3.5x / 4.5x / 5.5x	3.5x / 4.5x / 5.5x	3.5x / 5.0x
Ocular Weight	2.5x: 28.8g 3.0x: 34.4g 3.5x: 36.4g	3.5x: 44.6g 4.5x: 50.4g 5.5x: 55.8g	3.5x: 32.8g 4.5x: 35.7g 5.5x: 35.2g	3.5x: 20.4g 5.0x: 22.4g
Housing Material	Aircraft-grade Aluminum	Aircraft-grade Aluminum	Aluminum, Plastic	Aluminum, Plastic
Ocular Lens Material	Optical glass with anti-scratch & anti-reflective coating			


# FRAME SPECIFICATIONS

\*Frames weighed without lens or oculars

Enova Brand	Airon™ Standard	Airon™ Small	Airon™ Radius	Airon™ Tortoise	Airon™ Sport
Image					
Frame Color					
Frame Weight	17.9g	15.1g	15.8g	20.8g	16.7g
Frame Material	Aluminum	Aluminum	Aluminum	Acetate, Aluminum	Aluminum
Carrier Rx Material	High Index Trivex® with anti-reflective coating				Nylon

## EXPLORE MORE !

Take your loupes to the next level with Enova's selection of accessories! Our magnetic lights, 4K cameras, and more are sure to bring a cutting edge feel to your practice. Let one of our representatives help equip you for success.



To learn more about how Enova can help you see better, so you can perform better, visit our website or contact us.

### Enova Illumination

O: 651.236.8857  
E: sales@goenova.com W: goenova.com  
A: 6305 Sandburg Rd, STE 100, Minneapolis, MN 55427



F20546 Rev F

Copyright © 2026 Enova Illumination, Inc. All rights reserved. All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.