

GALILEAN

VIZIX® LOUPES



Prescription-friendly telescopes and carrier lens to create a precision fit unique to your eyes

Every frame comes with a magnetic adapter to securely hold a light or camera

All-glass telescopes lenses for high quality optical resolution

Built with all aluminum housings for easy cleaning and greater durability than plastic loupes

CUSTOM BUILT FOR YOU:

Bring precision to your profession with lightweight, aircraft-grade aluminum loupes from Enova's **Vizix® Galilean Loupe** collection. Performance and personalization go hand in hand with custom measured working distances, tailored declination angles, and optional lens prescriptions – all of which can be further complemented by an array of loupe and frame colors giving you a distinct look! Choose from Enova's Airon® family of frames or modern and Oakley® select frames.

MAGNIFICATION 2.5x, 3.0x, 3.5x

Select the level of magnification that best fits your needs. Our experienced representatives can help you identify a magnification level that is right for you.

CUSTOM OCULAR COLORS

Find the color combination that suits your style and personality! Available in silver, green, black, purple, and blue.

FUNCTIONAL FRAMES

Enova frames provide the complete package: modern style and superior durability. Comfort, ergonomics, and functionality are all seamlessly integrated into an attractive design.

OCULAR SPECIFICATIONS

*Oculars weighed as a pair without frame

	Galilean	Prismatic	Ergo Elite	Micro Ergo Elite
Image				
(Ocular Color)	 *3.5x is only available in Black/Silver			
Magnification	2.5x / 3.0x / 3.5x	3.5x / 4.5x / 5.5x	3.5x / 4.5x / 5.5x	3.5x / 5.0x
Ocular Weight	2.5x: 28.8g 3.0x: 34.4g 3.5x: 36.4g	3.5x: 44.6g 4.5x: 50.4g 5.5x: 55.8g	3.5x: 32.8g 4.5x: 35.7g 5.5x: 35.2g	3.5x: 20.4g 5.0x: 22.4g
Housing Material	Aircraft-grade Aluminum	Aircraft-grade Aluminum	Aluminum, Plastic	Aluminum, Plastic
Ocular Lens Material	Optical glass with anti-scratch & anti-reflective coating			

FRAME SPECIFICATIONS


*Frames weighted without lens or oculars

Enova Brand	Airon™ Standard	Airon™ Small	Airon™ Radius	Airon™ Tortoise	Airon™ Sport
Image					
Frame Color					
Frame Weight	17.9g	15.1g	15.8g	20.8g	16.7g
Frame Material	Aluminum	Aluminum	Aluminum	Acetate, Aluminum	Aluminum
Carrier Rx Material	High Index Trivex® with anti-reflective coating				Nylon

Oakley® Brand	Flak® 2.0 XL	CenterBoard
Image		
Frame Color		
Frame Weight	16.7g	12.9g
Frame Material	O-Matter, Unobtanium	Plastic
Carrier Rx Material	High Index Trivex® with anti-reflective coating	

EXPLORE MORE !

Take your loupes to the next level with Enova's selection of accessories! Our magnetic lights, 4K cameras, and more are sure to bring a cutting edge feel to your practice. Let one of our representatives help equip you for success.



To learn more about how Enova can help you see better, so you can perform better, visit our website or contact us.

Enova Illumination

O: 651.236.8857
E: sales@goenova.com W: goenova.com
A: 6305 Sandburg Rd, STE 100, Minneapolis, MN 55427



F20619 Rev E

Copyright © 2026 Enova Illumination, Inc. All rights reserved. All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.