

DENTAL LIGHT

AXIS ELITE / ATOM PLUS



ADVANCED COMFORT AND VISUALIZATION:

Enova developed the **Axis™ Elite LED Light** to be a lightweight, ergonomic illumination solution compatible with all Enova loupes.

The Axis Elite is one of the lightest dental lights on the market and the pivot connection allows for personalized alignment with field of view. The light can be paired with the Activ Elite neck battery or PP1.5 Elite hip battery to provide variable light intensity, and maximize magnification.

When used in conjunction with Enova's lightweight Vizix® loupes, clinicians benefit from superior illumination, visualization, and comfort, allowing them to perform with ease and confidence.

A TETHER-FREE LIGHT SOLUTION:

Experience the freedom of a wireless, tangle-free work environment with the **Atom Plus Wireless LED Light**. This lightweight, versatile solution provides a bright, white spotlight.

With a magnetic loupe adapter, it quickly and easily fits your existing loupes or most major brands with the perfect accessory. The reliable magnet attachment will provide all day hold. The Atom Plus wireless light is ideal for professionals who need a powerful lighting option for clinical or minor surgical procedures.

The package comes with two rechargeable batteries, a charger, detachable filter, and screwdriver- all packed in a carrying case.




TECHNICAL DATA

	Axis Elite	Atom Plus
Image & Color Options	 Not to size	 Not to size
Light Source	High Performance LED	LED Refractive Technology
Brightness	42,000 Lux @ 10 in	35,000 Lux @ 10 in
Color Temp.	5,700K	6,000K
Light Spot Size	3.25 in @ 12 in	4 in @ 12 in
Housing Material	Plastic	Aluminum
Compatible Battery	Activ Elite (Neck Battery) PP1.5 Elite (Hip Battery) <small>Not interchangeable</small>	Wireless Magnetic Battery
Battery Run-Time	Activ Elite: 5+ hours / PP1.5 Elite: 10+ hours	3+ hours
Light Weight	0.28oz, 7.8g (w/o filter)	0.64oz, 18.0g (w/o filter)
Battery Weight	Activ Elite: 5.3oz (150g) / PP1.5 Elite: 6.5oz (185g)	0.8oz (23g)
Battery Size	Activ Elite: 1x1 in, 21 in long / PP1.5 Elite: 4x2.5x1 in	2.2x1.1x0.6 in

	Axis Elite	Atom Plus
FEATURES		
Compatible adapters to fit most manufacturers' loupes	●	●
UV filter (Non-detachable) (Removing UV filter is optional)	●	●
Touch control - Brightness (Not available with PP1.5 Elite batteries)	●	●
Premium Package (Selected battery x2)	●	●
Standard Package (Selected battery x1)	●	-
WARRANTY		
LED covering manufacturer defects	3-year	3-year
Batteries and accessories	1-year	1-year

EXPLORE MORE !

Taking your light to the next level with Enova loupes is sure to bring a cutting edge feel to your practice. Let one of our representatives help equip you for success.



Connect with your sales representative now to experience a whole new light.

Enova Illumination

O: 651.236.8857
 E: sales@goenova.com W: goenova.com
 A: 6305 Sandburg Rd, STE 100, Minneapolis, MN 55427



F20676 Rev A

Copyright © 2026 Enova Illumination, Inc. All rights reserved. All trademarks are the property of Enova Illumination, Inc. unless otherwise noted.